# **NEWAR 2.0**

#### **Usage/Overview Guide**

#### August 2022

Jake Lorkin (jake@lorkin.cc)



#### Lab Room Setup





#### **POWER SWITCH**

Power board for computer and robot, ensure red light is on. Check room E-Stop if not.

### Robot Setup

• Close cage door, check door limit switch is properly closed.





• Check all 3 motor limit switches are in this position.

• The switch should 'click' when the motor moves to its end positions.

• Bump in motor coupler should be below limit switch roller.

## Robot Setup

• Check nothing is obstructing robots' range of movement, and cage area is clear.



### Power On

- Turn on motors by flicking motor power circuit breaker.
- Motors will energise and green LED's should be visible on each motor.
- Otherwise check control box is plugged in and has power from the wall.





### Power On

• Power on Control Box switch. Check control box e-stop is open.





#### **Control E-Stop**

Disables motor movement, motors will also be disabled if cage door is opened.



- Make sure control box is plugged into computer's USB port.
- Teknic ClearPath SC Hub should be visible on device manager.
- If not check all power connections and switches.

• Ensure door switch is closed and E-Stop open.

### Start Program-Basic

- Log into lab PC with staff or student account.
- All NEWAR files are saved locally under

#### 🚽 🗧 🛛 NEWAR Parallel Robot ~ File Home Share View ← C:\Users\Public\NEWAR Parallel Robot Name 📌 Quick access NUWAR\_Master\_Old 🔜 Desktop Software 🚽 Downloads \* ac\_clearpath-mc-sd\_manual 👸 Documents \* ClearPath\_MSP\_2.0.119 📰 Pictures \* Clearpath-SC User Manual MAINS MOTOR WIRING 👌 Music \* motor\_setup NEWAR Code NEWAR 2 Guide NEWAR Parallel Rok

#### C:\Users\Public\NEWAR Parallel Robot

• Run MATLAB and using the "Open" button open "NEWar.m"



• This line opens a '.coords' file that contains XYZ, coordinates and speed/acceleration info for the robot to follow

```
fid = fopen('demo.coords');
textLine = fgets(fid); % Read first line.
```

### Start Program - Basic

• A '.coords' file is a csv file that follows the structure

(X coord), (Y coord), (Z coord), (accel)\*, (vel)\*, (delay)\*

- \* means optional and all speed values are in RPM/s, delay is in milliseconds and coords in meters with the steel pole at the top being the origin.
- Delay is the minimum time a move will take, so once a move is done the remaining time will be waiting.
- optional parameters carry over onto previous lines until overwritten.
- image: constraint of the state of the s
- Save your coordinate file and run the MATLAB program.



• Motors should home, make a few clicking noises and then start the program.

Max Values:

- Motor Acceleration: 2000 RPM/s/s
- Motor Velocity: 500 RPM/s
- X: (0.7, -0.7). Y: (0.7, -0.7). Z: (-0.4, -0.7)